





LOOPS

Asas / Inoculateurs / Anse



HARD INOCULATION LOOP
asa de inoculación tipo rígido / inoculateur type rigide / ansa per inoculazione tipo rigido


STERILE
estéril / stérile / sterile





 **POLYSTYRENE**
poliestireno / polystyrène / polistirolo







SOFT INOCULATION LOOP
asa de inoculación tipo flexible / inoculateur type flexible / ansa per inoculazione tipo morbido

STERILE
estéril / stérile / sterile

 **POLYSTYRENE**
poliestireno / polystyrène / polistirolo

	CODE	CAPACITY	PACKAGING
	30010	1 µl (dark green)	BAG OF 20 PIECES bolsa de 20 piezas / sachet de 20 pièces / sacchetto da 20 pezzi
	30020	10 µl (blue)	BAG OF 20 PIECES bolsa de 20 piezas / sachet de 20 pièces / sacchetto da 20 pezzi
	30013	1 µl (dark green)	BAG OF 10 PIECES bolsa de 10 piezas / sachet de 10 pièces / sacchetto da 10 pezzi
	30023	10 µl (blue)	BAG OF 10 PIECES bolsa de 10 piezas / sachet de 10 pièces / sacchetto da 10 pezzi

	CODE	CAPACITY	PACKAGING
	30015	1 µl (light green)	BAG OF 20 PIECES bolsa de 20 piezas / sachet de 20 pièces / sacchetto da 20 pezzi
	30025	10 µl (light blue)	BAG OF 20 PIECES bolsa de 20 piezas / sachet de 20 pièces / sacchetto da 20 pezzi
	30016	1 µl (light green)	BAG OF 10 PIECES bolsa de 10 piezas / sachet de 10 pièces / sacchetto da 10 pezzi
	30026	10 µl (light blue)	BAG OF 10 PIECES bolsa de 10 piezas / sachet de 10 pièces / sacchetto da 10 pezzi